Nord GEL

EN	Phosphoric acid etching gel
DE	Ätzel mit Phosphorsäure
BG	Гел за ецване на фосфорна киселина
FR	Gel de mordantage à l'acide phosphorique
EL	Τζελ χάραξης φωσφορικού οξέος
IT	Gel mordenzante con acido fosforico
PL	Zel do wytrawiania kwasem fosforowym
ES	Gel de grabado de ácido fosfórico
TR	Fosforik asit aşındırma jeli

INSTRUCTION FOR USE EN

DESCRIPTION

Nord GEL is 37% thixotropic phosphoric acid gel with optimum level of viscosity. Syringe presentation with special tips allows placement of the etching gel in any desired location.

COMPOSITION

Diluent 40-60%, phosphoric acid 37%, silicon dioxide 10-20%, thickener 5-10%, pigment <1%.

Nord GEL does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

pH	<1
----	----

INTENDED PURPOSE AND CLINICAL BENEFITS

Nord GEL is intended to be used for tooth surface preparation procedures.

CLINICAL INDICATIONS

- **Tooth enamel and / or dentin etching.**

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, Nord GEL may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body - mouth. Tissues or body fluids contacted by the device - tooth.

INTENDED USER

Nord GEL is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental etchants. There is no need for specific training.

STERILITY

Nord GEL is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

Nord GEL is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of etchant is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as application tips, are supplied with the device.

INSTRUCTION FOR USE

PREPARATION:

1. Prepare surface as always. The working area must be kept dry.
2. If the surface to be etched is not freshly exposed, cleaning with oil-free prophylaxis paste such as Nord FASTE is recommended.
3. Attach application tip to the syringe.
4. Flow small amount of gel on mixing pad.

ETCHING:

1. Delicate push on plunger and apply layer of etch on all surface to be etched. Do not force plunger.
2. Leave etch gel in place for 15 seconds (dentine), 30 seconds (enamel).
3. Rinse with water and dry with air.
4. Verify etched surface is clean and dry.
5. Use an adhesive such as Nord BONDING LC.

WARNINGS

After the desired amount of material extruded, immediately remove application tip and close the syringe cap. Do not use Nord GEL for patients who have a history of

severe allergic or irritation reactions to product or to any of the ingredients. Nord GEL does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Subsequently, have the casualty drink 1 - 2 glasses of water immediately. Do not try to neutralize with bases, no charcoal. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If noticed any changes in the performance of the device, immediately rinse applied gel with copious amount of water, use suction system and do not let enter the mouth and swallow. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center.

SHELF-LIFE

Shelf-life of Nord GEL is 3 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Nord GEL is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

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